

Tutorial Four: Increased Affinity to a Peptide using Relax Design, Sequence Logo


Estimated time: 30 minutes

This tutorial will instruct users how to use Relax Design to increase the affinity between two peptides, including how to assess whether affinity has increased. The Sequence Logo tool will also be introduced.

In order to run this tutorial, you will need to know how to use the tools from Tutorial One, Two, and Three.

(Caveat: The estimated time does not include the time it takes for jobs to run on the server. There will be several pauses while you wait 20 minutes or more for structures to run on the cloud before you can move on to more steps.)

1) Run Relax & Minimize



For this tutorial the structure 10J5 will be used to model how Design can be used to increase peptide affinity. First the 10J5 structure will undergo Relax and Minimize to find a starting structure with optimized energy.


Select a Project or Create a Project where you will run this tutorial. Load structure 10J5. Run Relax, each round with 5 repeats, until the lowest energy conformation has been reached.

Minimize this final energy structure and rename the file *MinRelaxed 1OJ5*. This process of Relax rounds and Minimize is described in detail in Tutorial Two.

2) Design Increased Affinity Between Protein and Peptide

To increase the affinity between the two peptide chains in the 1OJ5 structure, residues involved in protein-protein interactions will be mutated. Using Relax Design, alternate amino acids will be sampled at selected residues to increase the affinity between the two peptide chains.

Click on *MinRelaxed 1OJ5* and click the  in order to visualize the protein in the Structure Viewer. One way to zoom in on the region of interest is to switch the selection mode from **Residues**, which is the default, to **Chains**. Then click on the small peptide and click  in order to center it. While the peptide is still selected, click  in order to show sticks as demonstrated below.



Switch the selection mode back to **Residues** and select the residues whose side chains are oriented towards the protein-peptide interface. It is important to note that prolines should not be selected for mutations because it is usually disruptive to mutate away from prolines. Click

 in order to show the sequences and highlight residues 106, 113, 116, and 117 (hold down command or control while selecting).

Click  and change the color in order to make it easier to see that these residues are important to binding. Create a Selector for the peptide selection by clicking , then clicking  to change the name to Peptide. Click the  to the left of the structure name in order to select it.


To sample alternate amino acids for the selected residues in the peptide click  in the right window. Choose  and enter . Then highlight all possible mutations by clicking  in the top left corner of the mutation table as shown below.


Design mutations

Pos.	Original Residues	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
106	L	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
113	L	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
116	L	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
117	L	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V

Click . Once the Design has run change the name of the Design run to *relD of 1OJ5* and click on the folder in order to bring the results into the center window.

Click  in the lower center window and change the button to  which will display which mutations were favorable. Click the  on the top left of the Sequence Logo window and click the arrow to the right of it to select the Peptide selector. This will hide all the residues that were not part of the selector. Scroll to the right to see the four residues in the selector as shown below.



The original sequence for these four positions was LLLL. You can see that L was only favorable at position 113. Position 106 preferred Q in all of the simulations. Positions 116 and 117 favored GA and I respectively. Results above are from a larger number of repeats in order to show diversity.

To look at the structure to see why these mutations were preferred switch

Sequence Logo ▲

back to

Structure Viewer ▲

. Now rank the structures by

Score

and click only the



for the top five structures.

Click in order to show their sequences. Then click in order to highlight the sequence differences. Scroll over to the right to see them.




In this example there were only two positions, 113 and 116 as highlighted in the sequences above, that differed across the top five structures. As suggested in the sequence logo, residue 113 would be most stable as an isoleucine or leucine and residue 116 as alanine or glycine.

The other two new mutations which were introduced at (Q106 and I117) are present across the top five structures suggesting they are best mutations for those positions.

Viewing and Comparing Potential Mutations

The mutations which were generated with Relax Design can be viewed and visually compared to the starting residues. The lowest energy structure from the Relax Design in this example is structure 25. Deselect the  for the four structures and select  to show structure 25.

Click on the whole structure and click  in order to show sticks. Select residues 106 to 119 by showing the sequence via  and clicking on the residues. Then click  and change the color in order to make it easier to see that these residues. Shown below are the three positions that were mutated in the structure.



With the L106Q mutation this side chain points away from the structure so it may be favorable because of improvement in solvation.

For L116A, the smaller side chain allows the peptide to shift towards the protein's beta sheet (see image to the right). The other preferred mutation at this site, glycine, also had a



smaller side chain to possibly allow the peptides to shift closer. This may not have been sterically possible before.

For L117I, the isoleucine packs into the gap which may have improved packing compared to leucine did.

In order to create an alignment of a specific region between two structures like the image in purple and green on the previous page, click **Compare**. In the drop down window, select the relaxed (*relID of 1OJ5*) and minimized wild type structure (MinRelaxed 1OJ5).

Ensure that the mutant and wild type structures are two different colors and align everything except the peptide. Do this by selecting residues 1 to 105 for both proteins. Click **Superpose** to open a drop down and click **Selected backbone atoms**.

As long as you leave the peptide out of the selection, this will align the protein so you can see how the peptide shifted when the mutations were modeled. You may want to go into full screen mode by clicking  in order to visualize the differences better.

Using Relax and Minimize to Compare Mutations

In order to verify that the mutations are energetically favorable first Relax and Minimize the structure generated from Design as follows. Click the to the left of the structure name and run a Relax with 5 repeats. Then run a Minimize on the most energetically favorable structure from Relax. Compare the energy of the Relaxed and Minimized mutant with the Relaxed and Minimized wild type as shown below:

minRelaxed pepMutant									
minimize 40 results Selected: 0 out of 1 structure(s) Load... Compare  Download selected									
#	<input type="checkbox"/>	Structure	Score	RMSD	Lennard-Jon	Backbone To	Sidechain To	Solvation	Hydrogen Bo
1	<input type="checkbox"/>	structure 1	-209.42	0.00	-383.89	-20.13	99.10	226.59	-128.87

MinRelaxed 1OJ5 Loaded in viewer: 0 Remove									
#	<input type="checkbox"/>	Structure	Score	RMSD	Lennard-Jon	Backbone To	Sidechain To	Solvation	Hydrogen Bo
1	<input type="checkbox"/>	structure 1	-205.02	0.01	-379.71	-19.46	96.41	223.67	-125.30


In the example above, the energy improved by 4.4 REUs from the minimized 1OJ5 structure to the minimized Design result. This is a significant difference for only 3 point mutations. However without testing it is not clear that the mutations have improved affinity. The mutations with lower energy structures may have improved stability within the peptide which doesn't affect affinity.

Predicting Peptide Stability Versus Peptide Affinity

One way to computationally predict whether the mutations are favorable due to peptide stability is to model the peptide without the protein. You will need to create and download a structure for 1OJ5 with only the smaller peptide chain present.

You can manually delete selections of the protein in a text editor or you can use a visualization program such as Pymol. The manipulated file needs to be a .pdb file with only the peptide residues present.

To download the original PDB file click the 1OJ5 folder to bring it into the center window. Click the to the left of the structure name in order to select it. Then click  to open a drop down window that allows you to select pdb, fasta, or scores. Select **PDB files**. Once only the smaller peptide chain residues remain the PDB file should look like this:



ATOM	ID	Element	Chain	Residue	X	Y	Z	B-factor	Occupancy	AltLoc
ATOM	1677	N	LEU	B 106	9.780	13.464	20.414	1.00	0.00	N
ATOM	1678	CA	LEU	B 106	9.540	14.311	21.629	1.00	0.00	C
ATOM	1679	C	LEU	B 106	10.845	14.936	22.141	1.00	0.00	C
ATOM	1680	O	LEU	B 106	11.733	14.228	22.632	1.00	0.00	O
ATOM	1681	CB	LEU	B 106	8.882	13.448	22.719	1.00	0.00	C
ATOM	1682	CG	LEU	B 106	8.499	14.104	24.050	1.00	0.00	C
ATOM	1683	CD1	LEU	B 106	7.448	15.179	23.863	1.00	0.00	C
ATOM	1684	CD2	LEU	B 106	7.983	13.011	24.992	1.00	0.00	C
ATOM	1685	1H	LEU	B 106	8.914	13.071	20.104	1.00	0.00	H
ATOM	1686	2H	LEU	B 106	10.169	14.028	19.686	1.00	0.00	H
ATOM	1687	3H	LEU	B 106	10.417	12.727	20.641	1.00	0.00	H
ATOM	1688	HA	LEU	B 106	8.867	15.123	21.357	1.00	0.00	H
ATOM	1689	1HB	LEU	B 106	7.963	13.023	22.318	1.00	0.00	H
ATOM	1690	2HB	LEU	B 106	9.557	12.630	22.970	1.00	0.00	H
ATOM	1691	HG	LEU	B 106	9.375	14.588	24.483	1.00	0.00	H
ATOM	1692	1HD1	LEU	B 106	7.203	15.620	24.829	1.00	0.00	H
ATOM	1693	2HD1	LEU	B 106	7.833	15.953	23.199	1.00	0.00	H
ATOM	1694	3HD1	LEU	B 106	6.551	14.739	23.428	1.00	0.00	H
ATOM	1695	1HD2	LEU	B 106	7.705	13.456	25.948	1.00	0.00	H
ATOM	1696	2HD2	LEU	B 106	7.111	12.530	24.548	1.00	0.00	H
ATOM	1697	3HD2	LEU	B 106	8.765	12.269	25.151	1.00	0.00	H
ATOM	1698	N	PRO	B 107	11.004	16.249	22.013	1.00	0.00	N
ATOM	1699	CA	PRO	B 107	12.223	16.872	22.517	1.00	0.00	C
ATOM	1700	C	PRO	B 107	12.396	16.717	24.041	1.00	0.00	C

Load this new structure into CAD as a new project. Create the three point mutations by running a single **Design** run, noting that the position numbers have been renumbered. You do not need to create a Selector for this. You can select the mutations manually in the mutation chart of the Design tab. Create these mutations:

Pos.	Original Residues	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
1	L	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
2	P	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
3	P	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
4	T	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
5	E	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
6	Q	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
7	D	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
8	L	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
9	T	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
10	K	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
11	L	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
12	L	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
13	L	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V
14	E	A	R	N	D	C	E	Q	G	H	I	L	K	M	F	P	S	T	W	Y	V

Leave the number of repeats as 1 and click . Rename the folder Peptide Mutant.

The WT Peptide and this new Peptide Mutant should be Relaxed and Minimized before comparing their energies. Since the peptide is so small, running 5  repeats is more than enough for both the WT and the new mutant. Then run  on the one with the best energy for both the WT and the new mutant. Using these, click  to see the energy difference.

#	Structure	Score	RMSD	Lennard-Jon	Backbone To Sidechain	To Solvation	Hydrogen Bc	
1	structure 1	-16.23	0.00	-32.61	-3.15	11.39	22.80	-10.26

#	Structure	Score	RMSD	Lennard-Jon	Backbone To Sidechain	To Solvation	Hydrogen Bc	
1	structure 1	-15.76	0.00	-32.58	-3.18	9.91	22.32	-9.34

In this case, the energy improved by 0.5 REUs, suggesting that the mutants are more likely to affect affinity to the protein than inner stability.